









	Batch No.		Do not re-sterilise
	Expiry date		CE mark and identification number of the Notified Body
	Sterilised by Ethylene oxide		Medical device
	MR CONDITIONAL (According to ASTM F2503)		Single sterile barrier with internal protection
	Unique device identification		Read the instructions for use available at https://www.mikai.us/downloads/
	Store in a cool, dry place		

Note: This device may have the following non-harmonised symbols (abbreviations) in the description:

- ST: indicates the sterility status of the device;
- Tl: indicates that the device is made of titanium;
- Fxx: the letter F followed by a number indicates the length of the main thread present in the device;
- Lx: the letter L followed by a number indicates the main length of the device;
- xxPZ: preceded by a number, indicates in the description the number of devices in the package.